Separation Techniques – Partitioning of Indicators



Background: A two-phase solvent system is made up of a polar and a non-polar solvent. Based on the differential solubility of the indicator added to such a two-phase system, and the subsequent partitioning of the chemical in one of the phases, we can determine whether the indicator is polar or non-polar.

This principle of partitioning can be used to separate components of a mixture using Partition chromatography wherein the mobile phase could be a mixture of solvents differing in their polarities and the stationary phase is also some solvent bound to a matrix such as paper.

In this experiment, we will be able to determine the nature of different indicators depending on their affinity for a particular solvent in a two phase solvent system. We will be using diethyl ether and distilled water as the two solvents to form the two-phase solvent system. Methyl orange, Sudan red and Thymol blue will be used as the indicators.

Category	Particulars
Chemicals	 Acetone Diethyl ether Distilled water Pre-made Methyl orange stain solution Pre-made Sudan red (IV) stain solution Pre-made Thymol blue stain solution
Apparatus/Facilities	 Plastic droppers Micropipette (1000 µL) and tips Test tubes Test tube stand Tip discard beaker containing tap water

Requirements:

Procedure:



- 1. Add 2 mL of distilled water to a test tube.
- 2. Then to the same test tube, add 2 mL of diethyl ether.
- 3. The test tube now contains a two-phase solvent system with diethyl ether and water forming the upper and lower phases, respectively.
- 4. Add a drop of methyl orange to the two-phase solvent system and shake the tube well. Based on its polarity, methyl orange will dissolve and be partitioned into one of the two phases.
- 5. Repeat the procedure for Sudan red and Thymol blue stain solutions in fresh tubes for each indicator, and record your observations.

Points for Discussion:

- 1. Identify the polar and non-polar stains based on their differential solubility in the two-phase solvent system.
- 2. Which of the three indicators can be used for staining waxy tissue in a plant leaf cross-section for microscopic observation?